

Clean Bench, Bio Reactor Type (이동식 접종기형 무균대)

Lab Safety

Application

이동식 접종기 및 생물반응기형 무균대는 배지접종, 식물 조직배양, 배지 작업, 의리기 조립, 의약품 작업 등 청정작업에 사용됩니다

Features


- ▷ 콘트롤러는 조작 버튼 형태로 사용과 상태 확인이 용이
- ▷ 초기 작동 시 진동 및 소음을 최소화한 브로워 모터의 연결부를 특수 처리 마감
- ▷ 평기류의 무진동 작업 테이블 타입



Bio reacting (Inoculation of a seed germ)



Horizontal non-vibration type clean bench **DK-2000H**



Dimensions **DK-2000H

Specifications

Model No.	DK-2000H	
Type	Scale bioreactor moving incubation & Horizontal non-vibration type	
**Dimensions, w x d x h / mm (In)	2000 x 900 x 1250	
(Out)	2180 x 1200 x 2000	
Purity	Lower than CLASS 100	
Main filter	H14 HEPA filter (EN 1822-1 compliant) Efficiency: 99.995% of the particles < 0.3 μm (Integral type, 2000 x 1200 x 75mm)	
Pre-filter	Nylon material filter	
Controller	Hour meter, Wind velocity HI/LOW switch, UV Lamp & FL Lamp S/W, Power S/W	
*Winding Velocity	0.3 ~ 0.5 m/sec	
Velocity Control	4 Stage (High, Middle, Low, Stop)	
Display	LCD Touch screen digital display (Speed, UV, Light, Used hour, Instruction manual)	
Material	(In)	Stainless steel plate (SUS 304)
	(Out)	Steel plate with powder coating finished (SCP-1)
	(Table, Non-vibration type)	(Top) Stainless steel (SUS304), 995 x 890 x 750mm x 2 sets / (Frame) Square steel pipe, 40 x 40mm with Foot master (40F)
	(Door)	Roll blind with Height adjustable wire / Sliding glass door, Tempered safety glass
Fluorescent lamp (FL)	40W x 4 pcs	
Germicidal lamp (UV)	40W x 2 pcs	
Utility	Consent box (220VAC/2 Gang), Gas cock, Vacuum cock, KPC One touch fitting/Air coupler, Hour meter (Built-in)	
Power source	220 VAC, 60 Hz, 1Phase	

*In accordance with SPS-KCLI2006-6235

주) 본 카탈로그의 사진, 내용 및 사양은 성능 및 품질개선을 위하여 예고 없이 변경될 수 있습니다.